



LoocidBCP™ Drill Protocols

LoocidBCP™ Drill Bits used for Nobel Biocare implants NobelActive

NB Active Ø	Loocid Starter Drill Bit (Optional)	Soft Bone IV/ Ø Loocid	Very Soft Bone/ Ø Active	Medium Bone/ Ø Loocid	Medium Bone/ Ø Active	Very Hard Bone/ Ø Loocid	Very Hard Bone/ Ø Active
NP 3.5mm	2.0mm	2.8mm	2.8mm	2.8mm	2.8 (3.2)mm	2.8mm	3.2mm
RP 4.3mm	2.0mm	2.8mm	2.8 (3.2)mm	3.2mm	3.6mm	3.2mm	3.6 (4.2)mm
RP 5.0mm	2.0mm	3.2mm	3.6mm	4.2mm	4.2mm	4.2mm	4.2 (4.6)mm
WP 5.5mm	2.0mm	3.2mm	3.6 (4.2)mm	4.2mm	4.2 (4.6)mm	4.2mm	4.6 (5.0)mm

The LoocidBCP™ Drill Bits can be used in conjunction with the Nobel Biocare Active Implant system. The surgical sequence is:

- 1) LoocidBCP™ Starter Drill Bit:
2.0mm (optional use only to verify drilling direction)
- 2) LoocidBCP™ Drill Bit:
The size of the LoocidBCP™ Drill Bit used after the optional Loocid 2.0mm Drill Bit or as the first drill bit for the implant bed preparation is shown in the table above
- 3) Nobel Biocare Active Drill Bit:
The Nobel Biocare Active Drill Bit following the use of the LoocidBCP™ Starter Drill Bit is shown in the table above. Please finish the implant bed preparation used for the correct implant diameter and bone density by following the regular Nobel Biocare Active Drill Bit protocol

