



# LoocidBCP™ Drill Protocols

## LoocidBCP™ Drill Bits used for Nobel Biocare implants NobelReplace CC and NobelReplace Tapered

NB Replace Ø	Loocid Starter Drill Bit (Optional)	Soft Bone/ Ø Loocid	Soft Bone/ Ø Replace	Medium Bone/ Ø Loocid	Medium Bone/ Ø Replace	Dense Bone/ Ø Loocid	Dense Bone/ Ø Replace
3.5mm	2.0mm	3.2mm	3.5mm	3.2mm	3.5mm	3.2mm	3.5 (3.5 Dense Bone Drill) mm
4.3mm	2.0mm	4.2mm	4.3mm	4.2mm	4.3mm	4.2mm	4.3 (4.3 Dense Bone Drill) mm
5.0mm	2.0mm	4.7mm	5.0mm	4.7mm	5.0mm	4.7mm	5.0 (5.0 Dense Bone Drill) mm
6.0mm	2.0mm	5.5mm	6.0mm	5.5mm	6.0mm	5.5mm	6.0 (6.0 Dense Bone Drill) mm

The LoocidBCP™ Drill Bits can be used in conjunction with the NobelReplace implant system. The surgical sequence is:

- 1) Optional Loocid Starter Drill Bit:  
2.0mm
- 2) LoocidBCP™ Drill Bit:  
The size of the LoocidBCP™ Drill Bit used after the optional Loocid 2.0mm Drill Bit or as the first drill bit for the implant bed preparation is shown in the table above
- 3) NobelReplace Drill Bit:  
The NobelReplace Drill Bit following the use of the LoocidBCP™ Starter Drill Bit is shown in the table above. Please finish the implant bed preparation used for the correct implant diameter and bone density by following the regular NobelReplace Drill Bit protocol

